

TN-1263

# Extraction and Analysis of Fentanyl and Analogs from Whole Blood Using a Kinetex™ F5 LC Column

Charissa Goggin<sup>1</sup> and Szabolcs Sofalvi<sup>2</sup>

<sup>1</sup>Phenomenex Inc., 411 Madrid Ave, Torrance, CA 90501 USA

<sup>2</sup>Cuyahoga County Medical Examiner's Office, Cuyahoga County Regional Forensic Science Laboratory, 11001 Cedar Avenue, Cleveland, OH 44106 USA



## Introduction

Fentanyl and fentanyl analogs are high potency synthetic opioids that are 80-100x more powerful than morphine. Although fentanyl use is strictly regulated, new fentanyl analogs keep emerging with molecular structures that are closely related, making separation, identification, and quantitation difficult for analysts. This method explores the identification of 26 compounds and quantitation of 24 compounds using multiple sample preparation techniques to clean-up the whole blood matrix. A core-shell LC column was utilized for its high efficiency and selective properties and the F5 phase was chosen for separating geometric and structural isomers.

## Methods

### Sample Pre-treatment

500 µL of a whole blood sample was lysed osmotically with 500 µL of DI Water and internal standard was added. The samples were vortexed for 7 seconds and they sat for 10 minutes. Protein precipitation was accomplished by drop-wise addition of 1.5 mL of ice cold methanol/acetonitrile (10:90) while vortexing. The samples were then incubated for 10 min in the freezer and then centrifuged at 3,500 rpm (2,800 xg) for 10 min. The supernatants were poured into new 16 x 125 mm test tubes followed by the addition of 3 mL DI water and vortexed for 7 sec.

## SPE Conditions

**Cartridge:** Strata™ Screen C (55 µm, 70 Å), 200 mg/6 mL

**Part No.:** 8B-S016-FCH

**Condition:** 3 mL of Methanol

**Equilibrate:** 2 mL of DI water

**Load:** Pre-treated samples were loaded onto the cartridge under 2 psi

**Wash 1:** 2 mL of DI water (All 3 wash steps were eluted under 5 psi positive pressure)

**Wash 2:** 1.5 mL of 1.0 M Acetic acid

**Wash 3:** 2 mL of Methanol

**Dry:** Cartridge for 10 minutes

**Elute:** 2x 1.5 mL of Dichloromethane/2-propanol/Ammonium hydroxide (78:20:2) under 5 psi. Full flow was applied at the end of each elution step for 30 sec.

**Dry Down:** Evaporated to dryness

**Reconstitute:** 40 µL of initial mobile phase. The samples were then vortexed for 7 sec and centrifuged at 2,800 xg for 5 min. The supernatants were transferred to conical inserts for UHPLC-MS/MS analysis.

## LC Conditions

**Column:** Kinetex 2.6 µm F5

**Dimensions:** 150 x 2.1 mm

**Part No.:** 00F-4723-AN

**Mobile Phase:** A: 0.1% Formic acid in Water  
B: 0.1% Formic acid in Acetonitrile

Gradient:	Time (min)	% B
	0	20
	17	25
	17.25	95
	18.75	95
	19.5	20
	21	20

**Flow Rate:** 0.5 mL/min

**Column Oven:** 40 °C

**Temperature:**

**Autosampler:** 5 °C

**Temperature:**

**Detector:** MS/MS (Thermo TSQ Quantis™ Triple Quad Tandem)


Analyte Name	Retention Time (min)
1. Norfentanyl	1.19
2. N-methyl norfentanyl	1.19
3. Norcarfentanil	1.47
4. Methoxyacetyl fentanyl	3.07
5. Acetyl fentanyl	3.51
6. β-hydroxy fentanyl	3.67
7. Benzyl fentanyl	3.95
8. para-methoxyacetyl fentanyl	4.67
9. Alfentanil	5.02
10. 4-ANPP	5.50
11. Acryl fentanyl	5.93
12. Fentanyl	6.42
13. para-fluoroacryl fentanyl	8.01
14. Cyclopropyl fentanyl	8.59
15. para-fluorofentanyl	8.82
16. 2-furanyl fentanyl	9.04
17. (E)-crotonyl fentanyl	9.35
18. (±)-trans-3-methylfentanyl	9.74
19. Fentanyl carbamate	10.27
20. (±)-cis-3-methylfentanyl	10.31
21. Carfentanil	11.40
22. Butyryl fentanyl	11.66
23. Sufentanil	13.76
24. para-fluoroisobutyryl fentanyl	14.29
25. Phenyl fentanyl	15.13
26. Cyclopentenyl fentanyl	15.59



## Results

Figure 1.

Separation and Analysis of Fentanyl and Metabolites Using Kinetex™ F5 LC Column



## Conclusion

Using whole blood as a matrix requires pre-treatment and sample preparation to clean-up prior to analysis. Employing a solid phase extraction method allows for matrix clean-up and the sample is concentrated which leads to more sensitive LC-MS/MS results. The Kinetex F5 LC column utilized a core-shell media to help which has 5 modes of selectivity: hydrophobic, aromatic, electrostatic, steric, and hydrogen bonding. This makes the phase extremely useful for isomeric analysis, and a great choice for fentanyl analog analysis.

## References

1. Szabolcs Sofalvi, Eric S Lavins, Ian T Brooker, Claire K Kaspar, John Kucmanic, Carrie D Mazzola, Christie L Mitchell-Mata, Cassandra L Clyde, Rindi N Rico, Luigino G Apollonio, Charissa Goggin, Brittany Marshall, Danielle Moore, Thomas P Gilson, (2019). Unique Structural/Stereo-Isomer and Isobar Analysis of Novel Fentanyl Analogues in Postmortem and DUID Whole Blood by UHPLC-MS-MS, Journal of Analytical Toxicology, Volume 43, Issue 9, bkz056, <https://doi.org/10.1093/jat/bkz056>

Phenomenex would like to thank Szabi Sofalvi and the rest of the Toxicology Department at Cuyahoga County Medical Examiner's Office for the use of their data and their time in developing this method.

## Ordering Information

## Kinetex™ LC Columns

5 µm Minibore Columns (mm)			SecurityGuard™ ULTRA Cartridges‡	
Phases	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
F5	<a href="#">00B-4724-AN</a>	<a href="#">00D-4724-AN</a>	<a href="#">00F-4724-AN</a>	<a href="#">AJ0-9322</a>
for 2.1 mm ID				

for 2.1 mm ID

5 µm MidBore™ Columns (mm)		SecurityGuard ULTRA Cartridges†	
Phases	100 x 3.0	150 x 3.0	3/pk
F5	<a href="#">00D-4724-Y0</a>	<a href="#">00F-4724-Y0</a>	<a href="#">AJ0-9321</a>

for 3.0 mm ID

5 µm Analytical Columns (mm)				SecurityGuard ULTRA Cartridges <sup>‡</sup>	
Phases	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	3/pk
F5	<a href="#">00B-4724-E0</a>	<a href="#">00D-4724-E0</a>	<a href="#">00F-4724-E0</a>	<a href="#">00G-4724-E0</a>	<a href="#">AJ0-9320</a>
					for 4.6 mm ID

for 4.6 mm ID

2.6 µm Minibore Columns (mm)				SecurityGuard ULTRA Cartridges <sup>‡</sup>	
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
F5	<a href="#">00A-4723-AN</a>	<a href="#">00B-4723-AN</a>	<a href="#">00D-4723-AN</a>	<a href="#">00F-4723-AN</a>	<a href="#">AJ0-9322</a>
					for 2.1 mm ID

for 2.1 mm ID

2.6 µm MidBore™ Columns (mm)			SecurityGuard ULTRA Cartridges†	
Phases	50 x 3.0	100 x 3.0	150 x 3.0	3/pk
F5	<a href="#">00B-4723-Y0</a>	<a href="#">00D-4723-Y0</a>	<a href="#">00F-4723-Y0</a>	<a href="#">AJ0-9321</a>
for 3.0 mm ID				

for 3.0 mm ID

2.6 µm Analytical Columns (mm)			SecurityGuard ULTRA Cartridges <sup>‡</sup>	
Phases	50 x 4.6	100 x 4.6	150 x 4.6	3/pk
F5	<a href="#">00B-4723-E0</a>	<a href="#">00D-4723-E0</a>	<a href="#">00F-4723-E0</a>	<a href="#">AJ0-9320</a>
for 4.6 mm ID				

for 4.6 mm ID

1.7 µm Minibore Columns (mm)			SecurityGuard ULTRA Cartridges†	
Phases	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
F5	<a href="#">00B-4722-AN</a>	<a href="#">00D-4722-AN</a>	<a href="#">00F-4722-AN</a>	<a href="#">AJ0-9322</a>
for 2.1 mm ID				

for 2.1 mm ID

† SecurityGuard ULTRA Cartridges require holder, Part No.: [AJ0-9000](#).

## Strata™ Screen C SPE

Format	Sorbent Mass	Part Number	Unit
Tube	100 mg	<a href="#">8B-S016-EAK</a>	1 mL (100/box)
	100 mg	<a href="#">8B-S016-EBJ</a>	3 mL (50/box)
	150 mg	<a href="#">8B-S016-SBJ</a>	3 mL (50/box)
	150 mg	<a href="#">8B-S016-SCH</a>	6 mL (30/box)
	200 mg	<a href="#">8B-S016-FBJ</a>	3 mL (50/box)
	200 mg	<a href="#">8B-S016-FCH</a>	6 mL (30/box)
	300 mg	<a href="#">8B-S016-RBJ</a>	3 mL (50/box)
	300 mg	<a href="#">8B-S016-RCH</a>	6 mL (30/box)
96-Well Plate	500 mg	<a href="#">8B-S016-HCH</a>	6 mL (30/box)
	50 mg	<a href="#">8E-S016-DGB</a>	2 Plates/Box



Have questions or want more details on implementing this method? We would love to help!  
Visit [www.phenomenex.com/chat](http://www.phenomenex.com/chat) to get in touch with one of our Technical Specialists

## Australia

t: +61 (0)2-9428-6444  
auinfo@phenomenex.com

## Austria

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

## Belgium

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

## Canada

t: +1 (800) 543-3681  
info@phenomenex.com

## China

t: +86 400-606-8099  
cninfo@phenomenex.com

## Czech Republic

t: +420 272 017 077  
cz-info@phenomenex.com

## Denmark

t: +45 4824 8048  
nordicinfo@phenomenex.com

## Finland

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

## France

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

## Germany

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

## Hong Kong

t: +852 6012 8162  
hkinfo@phenomenex.com

## India

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

## Indonesia

t: +62 21 3952 5747  
indoinfo@phenomenex.com

## Ireland

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

## Italy

t: +39 051 6327511  
italiainfo@phenomenex.com

## Japan

t: +81 (0) 120-149-262  
jpinfo@phenomenex.com

## Luxembourg

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

## Mexico

t: 01-800-844-5226  
tecnicomx@phenomenex.com

## The Netherlands

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

## New Zealand

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

## Norway

t: +47 810 02 005  
nordicinfo@phenomenex.com

## Poland

t: +48 22 51 02 180  
pl-info@phenomenex.com

## Portugal

t: +351 221 450 488  
ptinfo@phenomenex.com

## Singapore

t: 800-852-3944  
sginfo@phenomenex.com

## Slovakia

t: +420 272 017 077  
sk-info@phenomenex.com

## Spain

t: +34 91-413-8613  
espinfo@phenomenex.com

## Sweden

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

## Switzerland

t: +41 (0)61 692 20 20  
swissinfo@phenomenex.com

## Taiwan

t: +886 (0) 0801-49-1246  
twinfo@phenomenex.com

## Thailand

t: +66 (0) 2 566 0287  
thaiinfo@phenomenex.com

## United Kingdom

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

## USA

t: +1 (310) 212-0555  
info@phenomenex.com

## 🌐 All other countries/regions Corporate Office USA

t: +1 (310) 212-0555  
www.phenomenex.com/chat

## www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

# BE-HAPPY™ GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)

Strata, Kinetex, SecurityGuard, Be-Happy, and MidBore are trademarks of Phenomenex. TSQ Quantis is a trademark of Thermo Fisher Scientific Inc. *FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.*  
© 2025 Phenomenex, Inc. All rights reserved.

Have questions or want more details on implementing this method? We would love to help!  
Visit [www.phenomenex.com/chat](http://www.phenomenex.com/chat) to get in touch with one of our Technical Specialists

 **phenomenex™**