

TN-0160

Fast, Simple Urine Drug Extraction with In-well Room Temperature Hydrolysis and Enzyme Removal to Extend LC Column Life

Amanda C. McGee¹, Caleb R. Schlachter, PhD¹, L. Andrew Lee, PhD¹, Stephanie Marin, PhD², Babak Sheikhpoor², and Bryan Tackett, PhD²

¹Integrated Micro-Chromatography Systems, Inc., 110 Centrum Drive, Irmo, SC 29063 USA

²Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA



Liquid chromatography with tandem mass spectrometry (LC-MS/MS) is an essential tool for detecting the presence of drugs and their metabolites in urine samples. However, highly water-soluble glucuronide conjugates can affect the sensitivity and specificity of these tests. Many laboratories use hydrolysis to cleave Glucuronic Acid from the parent compound to overcome this challenge, thereby improving analyte detection by LC-MS/MS. One widely used hydrolysis technique in urine drug testing is enzymatic hydrolysis using β -glucuronidase enzymes.

Next-generation β -glucuronidase enzymes, such as IMCSzyme® RT, are highly active at room temperature and require less enzyme for complete hydrolysis. This reduction in enzyme quantity puts less protein on column and prolongs column life. However, due to endogenous metabolites in human urine samples, not all next-generation β -glucuronidase enzymes can hydrolyze glucuronidated metabolites equally. Previously, IMCS has reported that naturally occurring components of urine, such as urea, as well as metabolites from diet, such as vitamins and flavonoids, can have a negative impact on β -glucuronidase hydrolysis efficiency and not all β -glucuronidases are equally affected.

In this technical note, we evaluate the hydrolysis efficiency of IMCSzyme RT by 1) hydrolyzing 4 human urine samples at a range of enzyme amounts, and 2) hydrolyzing 96 human urine samples at a fixed amount of enzyme. We also evaluated how sample treatment using $\beta\text{-}Gone^{\text{\tiny TM}}$ Plus $\beta\text{-}glucuronidase$ removal plates can improve column life and peak quality of analytes compared to Dilute-and-Shoot treatment methods, using a Kinetex $^{\text{\tiny TM}}$ 2.6 μm Biphenyl LC column

Sample Preparation

Urine samples were a subset of those submitted for drug screening. Selected samples were screened and confirmed positive for Morphine, Oxymorphone, Hydromorphone, Codeine, and Temazepam. Certified Drug-Free Urine (DFU) was used as a control sample and fortified with glucoronidated drug standards. IMCSzyme RT was used for hydrolysis of the samples and protein concentration was based on A280 resulting in approximately 2 μ g/ μ L. Hydrolyzed samples were diluted to 40 % Methanol with 5 % Formic Acid in Methanol elution solvent and filtered through a β -Gone Plus plate (Part No.: 8E-S323-TGA) to remove residual hydrolysis enzyme. Samples were diluted to 10 % Methanol with Water prior to injection on LC-MS/MS.

Hydrolysis with a range of enzyme amounts: 12 human urine samples were pooled into 4 subsets to make samples positive for specific analytes: (1) Morphine, Hydromorphone, Buprenorphine, Norbuprenorphine, Temazepam, Oxazepam; (2) Oxymorphone and Hydromorphone; (3) Hydromorphone; and (4) Codeine. Control samples were fortified with the same glucuronidated metabolites at 100 ng/mL of free base. Control and human samples were hydrolyzed at room temperature (~20 °C) for 15 minutes with IMCSzyme RT at multiple enzyme concentrations ranging from 0 – 200 μg (Table 1).

(Gone


Hydrolysis with 20 μ g of enzyme: 96 human urine samples positive for at least one of the analytes (Morphine, Oxymorphone, Hydromorphone, Codeine, O-Desmethyl Tramadol, Tapentadol, Norbuprenorphine, Buprenorphine, Oxazepam, and/or Temazepam) were hydrolyzed with IMCSzyme RT using 100 μ L of urine, 20 μ g of enzyme and 300 μ L of room temperature buffer. 20 μ L of internal standard were also added to each sample.

IMCSzyme RT Dilute-and-Shoot Samples: Based on ≥ 80 % analyte recovery hydrolysis results 20 μg of IMCSzyme RT were used per reaction, mixed with 100 μL of urine, 300 μL of buffer and 20 μL of Methanol to represent internal standard. Samples were diluted to 500 μL final volume with 10% Methanol and 0.5% Formic Acid.

Table 1. IMCSzyme RT Master Mix Preparations.

Urine Sample (μL)	IMCSzyme RT at 2 μg/μL (μL)	IMCSzyme RT in Reaction (µg)	Water (μL)*	Room Temperature Hydrolysis Buffer (μL)	Internal Standard (μL)
	0	0	50		
	2	4	48		20
100	4	8	46		
	6	12	44	300	
	8	16	42	300	
	10	20	40		
	30	60	20		
	50	100	0		

^{*}Water was added after hydrolysis, so every sample had equivalent volume



MS/MS Conditions

Electrospray: 1000 V
Sheath Gas: 55 arb
Auxiliary Gas: 11 arb
Sweep Gas: 1 arb
Ion Transfer Tube Temperature: 300 °C
Vaporizer Temperature: 300 °C

2.35

Flow Rate: 0.5 mL/min

Injection Volume: 10 μL

Temperature: 40 °C

Gradient: Time (min)

0

0.5

1.5

1.55

1

Instrument: Thermo Scientific™ Vanquish™ HPLC

B: 0.1 % Formic Acid in Acetonitrile

%В

5

5

95

95

5

5

Detection: MS/MS

Detector: Thermo Scientific TSQ Endura™ MS

Table 2. MRM Transitions

Table 2. With Transitions											
Analyte	Retention Time (min)	Precursor (m/z)	Product (m/z)	Collision Energy (V)	RF Lens (V)	Analyte	Retention Time (min)	Precursor (m/z)	Product (m/z)	Collision Energy (V)	RF Lens (V)
Oxymorphone- Glucuronide	1.10	478.183	284.040 460.097	28 22	142	Norbuprenorphine- Glucuronide	3.34	590.304	396.208 414.222	37 35	156
$Morphine3\beta\text{-}Glucuronide$	1.12	462.183	211.012 286.040	42 29	144	Tapentadol	3.79	222.183	107.111 121.111	26 21	105
Hydromorphone-3- Glucuronide	1.23	462.183	184.986 286.054	47 28	168	Tapentadol-D₃	3.79	225.183	107.111 396.222	26 55	102
Morphine	2.09	286.4	157.100 165.100 183.100	39 40 32	136	Buprenorphine- Glucuronide	4.07	644.243	414.222 468.262 340.110	49 40 30	281
Morphine-D ₃	2.09	289.1	165.054	41	119	Norbuprenorphine	4.15	414.27	396.208	27	162
Oxymorphone	2.29	302.152	227.058 284.058	29 19	108	Norbuprenorphine-D ₃	4.15	417.274	343.165 241.000	29 38	164
Oxymorphone-D ₃	2.29	305.152	287.040	20	106	Oxazepam-Glucuronide	4.54	463	287.000	13	111
Naloxone Glucuronide	2.34	504.17	310.169 328.222	28 21	96	Lorazepam Glucuronide	4.63	498.7	276.875 322.889	37 13	109
Hydromorphone	2.54	286.122	486.222 157.000 184.986	20 42 30	127	Amitriptyline-Glucuronide	4.81	454.183	233.040 278.111 255.000	19 22 36	130
Hydromorphone-D ₃	2.54	289.183	184.929	30	130	Temazepam-Glucuronide	4.83	477.152	301.058	13	115
Naltrexol Glucuronide	2.77	520.27	326.468 501.994	21 38	229	Buprenorphine	4.9	468.365	396.151 414.222	40 35	194
Dihydrocodeine- Glucuronide	2.83	478.183	199.040 302.111	47 31	151	Buprenorphine-D ₄ Oxazepam	4.9 5.34	472.274 287	400.222 241.000	40 23	186 104
Codeine-6-Glucuronide	2.85	476.213	282.169 300.111	26 30	160	Oxazepam-D ₅	5.34	292	269.000 245.986	15 23	106
Naloxone	2.91	328.22	267.879 310.022	24 17	126	Amitriptyline	5.39	278.243	117.111 233.111	23 17	96
Naloxone-D ₅	2.91	333.2	315.151	21	106	Amitriptyline-D ₃	5.39	281.183	233.040	18	103
Dihydrocodeine	2.93	302.183	199.058 201.058	33 29	122	Lorazepam	5.42	322.091	275.960 303.875	22 15	109
Dihydrocodeine-D ₆	2.93	308.239	202.000	34	139	Lorazepam-D ₄	5.42	326.7	280.946	23	106
Codeine	2.99	300.19	165.111 215.111	42 25	129	Temazepam	5.9	301.091	255.058 282.986	22 14	99
Codeine-D ₆	2.99	306.183	218.111	26	128	Temazepam-D ₅	5.9	306.091	260.040	23	98
Naltrexol	3.17	343.865	254.240 326.169	31 21	107	cTHC Glucuronide	6.29	521.335	327.111 345.111	22 12	128
Naltrexol-D ₃	3.17	347.24	329.240	21	128	cTHC	7.14	345.274	299.111	19	132
Tapentadol Glucuronide	3.24	398.183	107.169 222.183	45 22	129	cTHC-D ₃	7.14	348.312	327.040 330.111	15 15	121
						26 53		5 .3.51L			

Results and Discussion

Ninety-six urine samples were hydrolyzed with IMCSzyme RT using 20 μ g of enzyme in the reaction for Dilute-and-Shoot and β -Gone workflows. Samples were considered positive (+) with \geq 25 ng/mL of drug and negative (-) with < 25 ng/mL of drug (**Table 3**).

Control and pooled urine samples were analyzed. Hydrolysis results were converted into % hydrolysis due to the differing concentration of drug in control and patient samples. Hydrolysis was considered complete at $\geq 80\%$ hydrolysis. 20 μg of enzyme was found to be sufficient for adequate hydrolysis of targeted analytes.

IMCSzyme RT with Dilute-and-Shoot and IMCSzyme RT processed with β -Gone Plus were compared. The injection sequence involved injecting one blank and three neat standard injections (10-minute run time). This was followed by injecting 96 Dilute-and-Shoot or β -Gone Plus treated samples. The sample injections were accelerated to 2.35 minutes to reduce run time and mobile phase waste, but column wash and column equilibration times remained equivalent between methods. This sequence was repeated 20 times for IMCSzyme RT Dilute-and-Shoot samples and the IMCSzyme RT processed with β -


Table 3. Summary of Drug Screening Results on 96 Urine Samples. Samples are Considered Positive at >25 ng/mL and Negative <25 ng/mL.


Drug	Positives	Negatives
Temazepam	7	89
Morphine	17	79
Oxymorphone	36	60
Hydromorphone	51	45
Codeine	5	91

Gone Plus. Each sequence finished with a blank and three neat standard injections for a total of 2008 injections for IMCSzyme RT Dilute-and-Shoot and IMCSzyme RT processed with β -Gone Plus.

All sample treatment methods started with a brand new Kinetex $2.6\,\mu m$ Biphenyl, $50\,x\,4.6$ mm column from the same batch to reduce variability. Total lon chromatograms (TICs) of Morphine were plotted. The full width at half max and asymmetry of Morphine peaks were calculated. Full width at half max should be as low as possible. Low values indicate a nice tall sharp peak while high values would indicate a short wide peak. Peak symmetry should be around 1, values higher than 1 indicate asymmetrical fronting and values less than 1 indicate asymmetrical tailing. The TICs of Morphine show that IMCSzyme R1 Dilute-and-Shoot samples caused significant retention time shifting (**Figure 1**). Full width at half max started out at 0.03, increased to 0.04 by 200 injections and stayed at 0.04 up to 2000 injections. Peak symmetry started out at 1.2 and fluctuated between 0.8 and 1.2 through 2000 injections.

Figure 1. Example Total Ion Chromatogram for Morphine.





Conclusions

The combination of IMCSzyme® RT and the β -Gone™ Plus plate provided fast, in-well room temperature hydrolysis and efficient removal of residual β -glucuronidase enzyme in urine samples. The total sample preparation time using β -Gone Plus and IMCSzyme RT was 20 minutes. This is significantly shorter than using traditional SPE or manual protein precipitation. Enzyme removal using β -Gone Plus resulted in cleaner samples and extended column life of the Kinetex™ 2.6 μ m Biphenyl LC column compared to Dilute-and-Shoot analysis. This workflow hydrolyzes even the most challenging urine samples and reduces the risk of LC column clogging and costly instrument downtime.

Ordering Information

Kinetex Analytical Columns

2.6 μm Analytic	al Columns (mm)						SecurityGuard™ ULTRA Cartridges‡
Phases	30 x 4.6	50 x 4.6	75 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	3/pk
EVO C18	00A-4725-E0	00B-4725-E0	_	00D-4725-E0	00F-4725-E0	00G-4725-E0	AJ0-9296
PS C18	<u>00A-4780-E0</u>	00B-4780-E0		00D-4780-E0	00F-4780-E0	00G-4780-E0	AJ0-8949
Polar C18	00A-4759-E0	00B-4759-E0		00D-4759-E0	00F-4759-E0		AJ0-9530
Biphenyl		00B-4622-E0		00D-4622-E0	00F-4622-E0		AJ0-9207
XB-C18		00B-4496-E0	00C-4496-E0	00D-4496-E0	00F-4496-E0	_	AJ0-8768
C18	00A-4462-E0	00B-4462-E0	00C-4462-E0	00D-4462-E0	00F-4462-E0	_	AJ0-8768
C8	_	00B-4497-E0	00C-4497-E0	00D-4497-E0	00F-4497-E0	_	<u>AJ0-8770</u>
HILIC	_	00B-4461-E0	00C-4461-E0	00D-4461-E0	00F-4461-E0	_	AJ0-8772
Phenyl-Hexyl	_	00B-4495-E0	00C-4495-E0	00D-4495-E0	00F-4495-E0		<u>AJ0-8774</u>
F5	<u>00A-4723-E0</u>	<u>00B-4723-E0</u>	_	00D-4723-E0	00F-4723-E0	_	<u>AJ0-9320</u>

for 4.6 mm ID

[‡]SecurityGuard ULTRA Cartridges require holder, Part No.: <u>AJO-9000</u>

β-Gone β-Glucuronidase Removal Products

Part No.	Description	Unit
8B-S139-TAK	1 mL Tubes, Recombinant Enzyme	100/Box
8B-S322-DAK	1 mL Tubes, Non-Recombinant Enzyme	100/Box
8E-S139-TGA	96-Well Plate, Recombinant Enzyme	1/Box
8E-S322-DGA	96-Well Plate, Non-Recombinant Enzyme	1/Box
8E-S323-TGA	96-Well Plate Plus 30 mg/well, Recombinant/Non-Recombinant Enzyme	1/Box
8E-S323-UGA	96-Well Plate Plus 60 mg/well, Recombinant/Non-Recombinant Enzyme	1/Box
8N-S323-TUK	2 mL Centrifuge Tubes, Recombinant and Non-Recombinant Enzyme	100/Box

Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/Chat.

t: +61 (0)2-9428-6444 auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301 anfrage@phenomenex.com

Belaium

t: +32 (0)2 503 4015 (French) t: +32 (0)2 511 8666 (Dutch) beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681 info@phenomenex.com

t: +86 400-606-8099 cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077 cz-info@phenomenex.com

Denmark

t: +45 4824 8048 nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063 nordicinfo@phenomenex.com

France t: +33 (0)1 30 09 21 10 franceinfo@phenomenex.com

t: +49 (0)6021-58830-0 anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162 hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400 indiainfo@phenomenex.com

Indonesia

t: +62 21 5019 9707 indoinfo@phenomenex.com

t: +353 (0)1 247 5405 eireinfo@phenomenex.com

Italy t: +39 051 6327511 italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262 jpinfo@phenomenex.com

Luxembourg t: +31 (0)30-2418700 nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226 tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700 nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951 nzinfo@phenomenex.com

Norway t: +47 810 02 005 nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72 pl-info@phenomenex.com

Portugal t: +351 221 450 488 ptinfo@phenomenex.com

Singapore t: +65 6559 4364 sginfo@phenomenex.com

Slovakia t: +420 272 017 077 sk-info@phenomenex.com

Spain

t: +34 91-413-8613 espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950 nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20 swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246 twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287 thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367 ukinfo@phenomenex.com

t: +1 (310) 212-0555 www.phenomenex.com/chat

All other countries/regions Corporate Office USA

t: +1 (310) 212-0555 www.phenomenex.com/chat

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com



Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/phx-terms-and-conditions-of-sale.

Trademarks

β-Gone, Kinetex, SecurityGuard, and BE-HAPPY are trademarks of Phenomenex. Thermo Scientific, Vanquish, and TSQ Endural are trademarks of Thermo Fisher Scientific. IMCSzyme is a registered trademark of Integrated Micro-Chromatography Systems, Inc. Surine is a trademark of DYNA-TEK Industries, Inc.

Disclaimer

SecurityGuard is patented by Phenomenex, U.S. Patent No. 6.162,362

CAUTION: this patent only applies to the analytical-sized guard cartridge holder, and does not apply to SemiPrep, PREP, or ULTRA holders, or to any cartridges. Phenomenex is in no way affiliated with Thermo Fisher Scientific, Integrated Micro-Chromatography Systems, Inc. or DYNA-TEK Industries, Inc.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures. © 2023 Phenomenex, Inc. All rights reserved.

